

Mimitakara®



@goodHearing P1

User Guide



1. ACCESSORIES	P 1
2. FEATURES	P 2
3. SET UP	P 3
3.1) Wearing instructions	P 3
3.2) Recharging instruction	P 5
3.3) Switching between modes	P 7
4. PAIRING BLUETOOTH	P 9
5. SAFETY & MAINTENANCE	P 12
5.1) Ear bud maintenance	P 12
5.2) Battery maintenance	P 13
5.3) LED Statuses	P 14
6. ADDITIONAL INFORMATION	P 15
7. FAQ	P 16
8. SPECIFICATION	P 19
9. WARNING AND REMINDER	P 20

1. ACCESSORIES

- 1 Main unit with built-in ear buds (M size)
- 2 2 x Additional pair of L & S Size ear buds
- 3 USB cable



2. FEATURES



3. SET UP

3.1) Wearing instructions

- 1) Securely fit the red earphone into the right ear and the blue earphone into the left ear.
- 2) Turn the unit ON (When the unit is turned on the LED Indicator Light will flash blue.
- 3) The **goodHearing P1** has 11 adjustable volume settings, adjust the volume to suit your preferred specification. When the hearing aid is turned off and turned back on it will resume from its previous setting from when it was turned off.

CAUTION : LISTENING AT A HIGH VOLUMES MAY DAMAGE YOUR HEARING.

4) We suggest clipping the device to the outside shirt pocket or collar to avoid it from rubbing against objects in order to reduce the static sounds (Ensure the microphone is not obstructed for optimum sound quality).



3.2) Recharging Instruction

1) Charging the battery before use is recommended but not required.

goodHearing P1 comes pre-charged however initial charge will vary.

2) 3+ hrs for full charge gives up to 10 hours of usage time.



3) When plugged into a power source



Solid RED LED Indicator : Charging




Solid BLUE LED Indicator : Fully charged

3.3) Switching between modes

1) There are 3 pre-set modes on the **goodHearing P1** (Standard Mode & Noise Cancelling Mode).

- Standard Mode :
Mostly used indoor where the environment is not too noisy.
- Noise Cancelling Mode :
Used in louder environment such as outdoor public areas, restaurants.
- High-level Noise Cancelling Mode :
Providing superior noise cancelling in even the loudest environments, such as in transportation

2) To switch between the 2 modes click the  button (located in the center of the **goodHearing P1** device). A beep will be heard when switching between the 3 modes to indicate which mode it is currently set it.

- Mode 1 : A single beep will be heard
(Standard Mode)
- Mode 2 : 2 beeps will be heard
(Noise Cancelling Mode)
- Mode 3 : 3 beeps will be heard
(High-level Noise Cancelling Mode)

4. PAIRING BLUETOOTH

- 1) Have both devices within arm length of each other
- 2) Turn on **goodHearing P1**
- 3) Press and hold the M for 3 seconds
- 4) The LED Indicator will flash alternating BLUE and RED colors
- 5) Turn on the Bluetooth function on your mobile device
- 6) Scan for nearby devices
- 7) Pair with **goodHearing P1**
- 8) When pairing has been successful your phone will show that the status is connected with **goodHearing P1**

Paired devices will remember each other unless pairing information is deleted, so there is no need to pair again. **goodHearing P1** can pair and remember 2 bluetooth devices.

Repeat the steps outlined above to set-up your **goodHearing P1** with another Bluetooth device. Note that when a 3rd Bluetooth device is paired the oldest pairing Bluetooth device information will be deleted.

goodHearing P1 will connect with the latest connected device.

If you want to use another paired Bluetooth device, switch off the Bluetooth function of the current Bluetooth device, and then connect the Bluetooth device to the hearing device by selecting **goodHearing P1** on the Bluetooth device list and press connect.

When the Bluetooth device that is connected to the **goodHearing P1** is not playing any sort of media, the device will automatically switch back to being its primary function until the connected device starts playing media or receives a Should the Bluetooth be disconnected due to range or

any other reason, reconnect the **goodHearing P1** with the device that you wish to pair it with.

Please scan QR code to download the APP or search "goodHearing" on APP Store / Google Play



goodHearing



iOS




Android

5. SAFETY AND MAINTENANCE

5.1) Ear bud maintenance

- Before each use : Clean your ear thoroughly.
- After each use : Remove all remaining debris from the ear bud after each use. Ensure that the hole on the side of ear bud is also free from any debris.



Use the cleaning brush to clean any earwax adhered on the hearing aid, especially in and around the hearing aid and ear bud nipple.

5.2) Battery maintenance

- Do not consistently discharge the battery completely before recharging as this may cause depletion of the batteries life.
- We suggest recharging **goodHearing P1** regardless of its battery status before sleeping at night.
- Do not attempt to remove the battery as this may cause damage to the **goodHearing P1** device or may cause injuries.
- Ensure the device is fully charged before removing it from charge.

5.3) LED Indicator status

- Solid RED : Charging
- Flashing BLUE : ON (working normally)
- Flashing RED and BLUE :
goodHearing P1 is searching for Bluetooth devices to pair with.
- Flashing PURPLE :
Connected to a device and outputting media audio.

6. ADDITIONAL INFORMATION

- Keep the device out of reach for children.
- Do not attempt to open, fix, replace or remove any of the components.
- Avoid excessive volume levels and prolonged periods of use.
- Avoid using the device while you are sleeping.
- Keep the device in dry areas.
- Keep the device away from high temperature areas.

7. FAQ

Question : No Sound

Answer :

- Turn on the device
- Increase the volume
- Clean out the earphone caps
- Check if the earphone cable damages
- Recharge the device

Question : Charging problems

Answer :

- Check that the USB is connected correctly
- Ensure that the power supply is normal
(wall outlet switch is on)

Question : Can't power on

Answer :

Out of battery and needs to be charged.

Question : Is it compatible for both Android and iOS systems?

Answer :

Yes, is compatible for both Android and iOS systems.

Question : There is a lag when I change from listening to media on my phone and back to the hearing aid function is this normal?

Answers :

Yes, there is a slight lag as it depends on the operating system you are currently running on your phone. There may be up to 5 second delay.

Question : Can I pick up phone calls with the hearing aid?

Answers :

Yes, you can, with Android and iOS system your

phones will have a default tone that the hearing aid will play as a ringtone.

Question : There is still no sound on my SONY phone.

Answers :

Please turn off the following 3 sound settings.

1. Keyboard touch tone
2. Screen lock sound
3. Vibration on touch

Question : There is still no sound on my HTC phone.

Answers :

Please turn off the following 4 sound settings.

1. Keyboard touch tones
2. Touch sounds
3. Screen lock sounds
4. Pull-to-refresh sound

8. SPECIFICATION

Max. Output Sound Level Gain

45 ± 5 dB

Continuous Use

Up to 10 hours per charge

Operating Temperature

+10°C ~ +40°C

Operating Humidity

30~85% RH

Storage Temperature

-10°C ~ +60°C

Storage Humidity

30~85% RH

Dimension

45 x 33 x 14 mm

Weight

22.4 g

9. WARNING AND REMINDER

- Do not listen for excessive time periods.
- Do not use the device at high volumes.
- Do not attempt to fix the unit yourself, take it to an authorized dealer.
- Avoid liquids near the device.
- Avoid overcharging the device for excessive time periods.
- Avoid discharging the battery for long periods.
- Keep away from reach of children.
- Keep the device dry.
- Keep away from high temperature areas.
- Keep away from direct sunlight.
- Ensure the earphones are fitted securely in the ear to reduce sound leak.

