

Mimitakara® JAPAN

www.mimitakara.com

Oral Cleaner

Welcome to the Mimitakara family and congratulations on your purchase of the Mimitakara Oral Cleaner. By purchasing this Oral Cleaner we are sure you are aware of the benefits of keeping your oral products clean and sanitized.

Mimitakara's Oral Cleaner can keep an array of your oral products clean and sanitized, such as dentures, orthodontics, mouth guards and whitening trays is now easier than ever before. With a simple push of a button our Oral Cleaner will take care of the rest. The device is made hassle free and will automatically shut down after 20 min of cleaning and sanitizing.

UV-C Light is effective in destroying nucleic acids and breaking apart germ DNA. With their DNA broken, they can't function or reproduce and the organism dies. UV-C sanitizing has been an accepted practice since the mid-20th century. In fact, the 1903 Nobel Prize for Medicine was awarded to Niels Finsen for his use of UV against tuberculosis. It has since been used primarily in medical sanitation and sterile work facilities. UV-C light penetrates the bacteria's cell and disrupts the DNA killing the pathogens.

USING DENTURE CLEANER

Before use, unscrew the battery door compartment and insert 3 "AA" batteries

1. Push the small white clip to open the lid of your Oral Cleaner. Remove the cup container and curved tray from the inside compartment by lifting the cup container straight up and out.
2. Set the curved tray aside.
3. Fill the cup container with water to the indicator line.
4. Place the cup container back into the Oral Cleaner. When seated, push down firmly for a tight fit.
5. Place items to be cleaned directly into the cup container so they are submerged in the water. If you are cleaning smaller items, you can place them on the curved tray and then insert into the cup container filled with water.
6. Close the lid of the Oral Cleaner.
7. Choose the cleaning method: (1) UV-C ONLY / (2) UV-C + SONIC CLEANING

UV-C CLEANING ONLY :

- Press the "UV Only" button to start. The indicator light will flash blue to indicate it is working. It will switch off automatically after 20 minutes.

UV-C + SONIC CLEANING :

- Press the "UV + Sonic" button to start. The indicator light will flash red and blue to indicate that it is working. The sonic vibration will turn off automatically 5 minutes, but the UV-C bacteria cleaning function will continue for an extra 15 minutes.

8. Manual Switch-Off: To interrupt either cleaning cycles, simply press and hold the "UV Only" or "UV + Sonic" button or simply just open the lid.
9. After your Mimitakara Oral Cleaner. has completed its cleaning cycle, remove the items and discard the water from the cup. Rinse and dry the cup so it's ready for the next use.

CAUTION

- Do not place unit where it can be in danger of falling into sinks, toilets or bathtubs.
- Always make sure your Oral Cleaner's lid is closed before using. The unit shuts down when the lid is opened.
- Should the unit be dropped or is damaged, please do not attempt to fix the unit.
- Remove and rinse the cup container and curved tray after each use.
- Never place the unit base under running or immerse it under water for any reason. To clean, simply wipe with a moist cloth as needed and keep dry.
- Do not look directly at the UV Light.

TROUBLE SHOOTING

- Press hard enough on the blue mode buttons?
- Are the batteries in need of replacement?
- Is the lid closed properly?
- The UV-C bacteria cleaning light needs time to heat up and activate, it may take up to 2 - 3 min to activate upon first use.

SPECIFICATION

Model : Oral Cleaner

Power : 3 x AA Batteries

AC Adapter (no included)

Activation Duration : 20 min

Dimension : 4" x 4.5" x 4.3" Inches

Weight : 0.66 lb

Made In Taiwan

Designed in Japan

